



COVER ALL TUBING ENDS WITH POLY BAG DUST COVERS.

EACH FINISHED CTC-1322 WILL BE PACKED IN A 4 MIL. POLY BAG.

TUBING TOLERANCES: 84 IN. OR > ±2.00"
 BETWEEN 60 AND 84 IN. ±1.00"
 BETWEEN 12 AND 60 IN. ±0.50"
 12 IN. OR < ±0.25"

INSTALL CL-01009 CABLE TIES AT EACH CONNECTION.

BILL OF MATERIALS			REV DATE	DESCRIPTION
*NOT SHOWN IN DRAWING			A	NEW DEVELOPMENT.
PART#	Q.P.A.	DESCRIPTION	03/29/06	
S4111	1.00	1 LITER BAG		
BC-1001-I01	1.00	1 LITER DIE	B	COMPANY NAME CHANGE TO LONZA, FILM DESCRIPTION CORRECTED TO PLATINUM ULTRAPAK, NO CHANGE TO DESIGN.
PL-01048	-	PLATINUM ULTRAPAK INNER FILM	04/24/07	
PL-01049	-	PLATINUM ULTRAPAK OUTER BARRIER FILM		
TB-01062	2.00	BAG PORT WITH 1/4" BARB		
TB-01130	1.50 FT.	1/4" I.D. X 3/8" O.D. PVC TUBING		
TB-01254	1.50 FT.	1/4" I.D. X 3/8" O.D. C-FLEX NAD TUBING		
CL-01021	2.00	PINCH CLAMP		
CN-01065	1.00	MALE MPC WITH 1/4" BARB		
CP-01068	1.00	FEMALE MPC CAP		
CN-01061	1.00	FEMALE LUER WITH 1/4" BARB		
CN-01225	1.00	SWABABLE VALVE		
CP-01055	1.00	MALE COVER CAP		

Lonza Walkersville, Inc.
 8830 Biggs Ford Road
 Walkersville, MD 21793
 Phone 301.898.7025

APPROVALS DATE
 DFM JWP 04/24/07
 PROD.

CUSTOMER
 MARKETING
 Engineering

PRICE
 -

Scale DRAWING NO. REV.
 - CTC-1322 B