



COVER ALL TUBING ENDS WITH POLY BAG DUST COVERS.

EACH FINISHED CTC-1326 WILL BE PACKED IN A 4 MIL. POLY BAG.

TUBING TOLERANCES: 84 IN. OR > ±2.00"
 BETWEEN 60 AND 84 IN. ±1.00"
 BETWEEN 12 AND 60 IN. ±0.50"
 12 IN. OR < ±0.25"

INSTALL CL-01009 CABLE TIES AT EACH CONNECTION.
 TUBING LENGTHS 36 IN. AND OVER WILL BE ROLLED AND SECURED WITH A CL-01009 CABLE TIE.

BILL OF MATERIALS

*NOT SHOWN IN DRAWING

PART#	Q.P.A.	DESCRIPTION
S4256	1.00	50 LITER BAG WITH FITMENTS
BC-0050	1.00	50 LITER DIE
PL-01048	-	PLATINUM ULTRAPAK INNER FILM
PL-01049	-	PLATINUM ULTRAPAK OUTER BARRIER FILM
TB-01062	3.00	BAG PORT WITH 1/4" BARB
TB-01023	12.00 FT.	1/4" I.D. X 3/8" O.D. PCMGS TUBING
TB-01254	1.50 FT.	1/4" I.D. X 3/8" O.D. C-FLEX NAD TUBING
CL-01021	3.00	SMALL PINCH CLAMP
CN-01055	1.00	FEMALE MPC CONNECTOR WITH 1/4" BARB
CP-01020	1.00	MALE MPC PLUG
CN-01065	1.00	MALE MPC CONNECTOR WITH 1/4" BARB
CP-01068	1.00	FEMALE MPC CAP
CN-01061	1.00	FEMALE LUER WITH 1/4" BARB
CN-01225	1.00	SWABABLE VALVE
CP-01055	1.00	MALE COVER CAP

REV DATE	DESCRIPTION
A 03/23/06	NEW DEVELOPMENT.
B 04/24/07	COMPANY NAME CHANGE TO LONZA, FILM DESCRIPTION CORRECTED TO PLATINUM ULTRAPAK, NO CHANGE TO DESIGN.

Lonza Walkersville, Inc.
 8830 Biggs Ford Road
 Walkersville, MD 21793
 Phone 301.898.7025

CUSTOMER		APPROVALS	DATE
TITLE	DRM	JWP	04/24/07
50 LITER	PROD.		
PRICE	CUSTOMER		
-	MARKETING		
Scale	Engineering		
-			
DRAWING NO.	REV.		
CTC-1326	B		